

## Collection Instructions: Please read all instructions carefully

### Instructions for Men:

- Ø Clean the tip of the penis with the wrapped antiseptic towelette provided.
- Ø Pass a small amount of urine into the toilet.
- Ø Collect the specimen by passing the middle portion of urine into the sterile container provided until it is half full. Then empty the bladder into the toilet.

### Instructions for Women:

- Ø Wipe the entire outer vaginal area with the wrapped antiseptic towelette provided.
- Ø Pass a small amount of urine into the toilet.
- Ø Collect the specimen by passing the middle portion of urine into the sterile container provided until it is half full. Then empty the remaining urine in the bladder into the toilet.

## Containers:

