

Women's Health Swab Collection Guide: Vaginosis (BV)/Vaginitis, Sexually Transmitted Infection (STI)







Specimens submitted from sources not noted as acceptable for the collection kit cannot be tested. Up to 4 tests can be run from ONE Aptima® Multitest Swab.

Test Name	Order Code	Specimen Source	Collection Container	Transport Temp
CANDIDA SPECIES GROUP, CANDIDA GLABRATA, TRICHOMONAS VAGINALIS, NAAT	LAB10625	Vaginal	APTIMA® MULTITEST	Room Temp
CHLAMYDIA TRACHOMATIS / NEISSERIA GONORRHOEAE (CT/GC) RRNA	LAB9350	Cervical	THINPREP® VIAL	Room Temp
		Endocervical	APTIMA® UNISEX	Room Temp
		Vaginal	APTIMA® MULTITEST	Room Temp
		Urine (first catch)	STERILE CONTAINER	Refrigerate
CULTURE, GENITAL (CULTURE AND GRAM STAIN)	LAB8711	Cervix, Endocervix, Vaginal	eSwab® or BBL CultureSwab™	Room Temp
HERPES SIMPLEX VIRUS TYPES 1 AND 2, NAAT (ANOGENITAL)	LAB10598	Anogenital	APTIMA® MULTITEST	Room Temp
HERPES SIMPLEX VIRUS TYPES 1 AND 2, NAAT (NON-GENITAL SWAB)	LAB9370	Non-genital	UTM® VIRAL TRANSPORT	Room Temp
MYCOPLASMA GENITAL PROFILE, NAA (SWAB)	LAB9394	Endocervical, Vaginal	APTIMA® UNISEX	Room Temp
				Room Temp
MYCOPLASMA HOMINIS/UREAPLASMA CULTURE	LAB9637	Vaginal	UTM/UVT	Refrigerate
		Urine (first catch)	STERILE CONTAINER	Refrigerate
TRICHOMONAS VAGINALIS (TV), NAAT	LAB10600	Cervical	THINPREP® VIAL	Room Temp
		Vaginal	APTIMA® MULTITEST	Room Temp
		Endocervical	APTIMA® UNISEX	Room Temp
VAGINOSIS SCREEN (GRAM STAIN)	LAB8712	Cervix, Endocervix, Vaginal	eSwab® or BBL CultureSwab™	Room Temp
VAGINOSIS, BACTERIAL (BV), NAAT	LAB10624	Vaginal	APTIMA® MULTITEST	Room Temp
VAGINITIS PANEL (AFFIRM)	LAB8755	Vaginal	BD Affirm™ VPIII Swab	Refrigerate
VAGINITIS PANEL (BV, CANDIDA 5 SP., TV) NAA	LAB10638	Vaginal	APTIMA® MULTITEST	Room Temp
VAGINITIS PANEL PLUS (BV, CANDIDA 5 SP., TV, CT/GC) NAA	LAB10639	Vaginal	APTIMA® MULTITEST	Room Temp
VIRAL CULTURE, GENERAL	LAB9438	Varies	UTM® VIRAL TRANSPORT	Refrigerate

Candida species group includes: *Candida albicans*, *Candida parapsilosis*, *Candida dubliniensis*, *Candida tropicalis* **Mycoplasma genital profile includes:** *Mycoplasma genitalium*; *M. hominis*; *Ureaplasma spp* **BV includes:** *Lactobacillus gasseri*, *Lactobacillus crispatus*, *Lactobacillus jensenii*, *Gardenerella vaginalis*, *Atopobium vaginae* **BD Affirm includes:** *Candida species*, *Gardenerella vaginalis*, *Trichomonas vaginalis* **Viral culture may detect:** Adenovirus, CMV, enteroviruses, HSV, influenza, mumps, parainfluenza, RSV, V-Z cell culture

Women's Health Swab Collection Guide: Vaginosis (BV)/Vaginitis, Sexually Transmitted Infection (STI)

Specimens submitted from sources not noted as acceptable for the collection kit cannot be tested. Up to 4 tests can be run from ONE Aptima® Multitest Swab.

APTIMA® MULTITEST	APTIMA® UNISEX	APTIMA® URINE	BD Affirm™ VPIII Swab	UTM® Viral Transport	BBL CultureSwab™
					
STERILE CONTAINER	THINPREP® VIAL	eSwab® (white cap) Regular Flocked Swab	eSwab® (purple cap) Regular Flocked Swab	eSwab® (blue cap) Flex Minitip Flocked Swab	eSwab® (green cap) Minitip Flocked Swab
