

## Women's Health SWAB Collection Guide Container/Source

| COLLECTION CONTAINER  |            | APTIMA®<br>MULTITEST<br>Swab<br>(orange label) | APTIMA®<br>UNISEX<br>Swab<br>(white label) | eSwab™<br>Transport<br>Swab         | BD Affirm™<br>VPIII Swab | Universal<br>Transport<br>Media<br>(UTM/UVT)      |
|---|------------|--|--|-------------------------------------|--------------------------|---|
| SPECIMEN SOURCE   |            | Vaginal  | Endocervical                               | Cervical,<br>Endocervix,<br>Vaginal | Vaginal                  | Genital,<br>Non-Genital,<br>Rectal,<br>Throat, NP |
| Test Name   | Order Code |  |  |                                     |                          |   |
| Candida Species Group*, Candida Glabrata, Trichomonas Vaginalis, NAAT | LAB10625   | ■  |  |                                     |                          |   |
| Chlamydia Trachomatis / Neisseria Gonorrhoeae (CT/GC) RRNA            | LAB9350    | ■  | ■  |                                     |                          |   |
| Herpes Simplex Virus Types 1 And 2, NAAT (Anogenital Lesion)          | LAB10598   | ■  |  |                                     |                          |   |
| Herpes Simplex Virus Types 1 And 2, NAAT (Non-Genital Swab)           | LAB9370    |  |  |                                     |                          | ■   |
| Mycoplasma Genital Profile, NAAT† (Swab)                              | LAB9394    |  | ■  |                                     |                          |   |
| Mycoplasma hominis/Ureaplasma Culture                                 | LAB9637    |  |  |                                     |                          | ■   |
| Trichomonas Vaginalis, NAAT   | LAB10600   | ■  | ■  |                                     |                          |   |
| Vaginosis Screen (Gram Stain)   | LAB8712    |  |  | ■                                   |                          |   |
| Vaginosis, Bacterial (BV), NAAT‡                                      | LAB10624   | ■  | ■  |                                     |                          |   |
| Vaginitis Panel (Affirm)§   | LAB8755    |  |  |                                     | ■                        |   |
| Vaginitis Panel (BV, Candida 5 Sp., TV), NAAT                         | LAB10638   | ■  |  |                                     |                          |   |
| Vaginitis Panel Plus (BV, Candida 5 Sp., TV, CT/GC), NAAT¶            | LAB10639   | ■  |  |                                     |                          |   |
| Viral Culture, General¥   | LAB9438    |  |  |                                     |                          | ■   |

\* *Candida species group* Includes *Candida albican*, *Candida parapsilosis*, *Candida dubliniensis*, *Candida tropicalis*  
 † *Mycoplasma genital profile* Includes *Mycoplasma genitalium*; *M. hominis*; *Ureaplasma spp.*  
 ‡ *Vaginosis Bacterial (BV)* Includes *Lactobacillus gasseri*, *Lactobacillus crispatus*, *Lactobacillus jensenii*, *Gardenerella vaginalis*, *Atopobium vaginae*  
 § *Vaginitis Panel BD Affirm™* includes *Candida species*, *Gardenerella vaginalis*, *Trichomonas vaginalis*  
 || *Vaginitis Panel* includes *Bacterial Vaginosis*, *Candida species group*, *Trichomonas vaginalis*  
 ¶ *Vaginitis Panel Plus* includes *Bacterial Vaginosis*, *Candida species group*, *Trichomonas vaginalis*, *Chlamydia trachomatis/Neisseria gonorrhoeae*  
 ¥ *Viral Culture, General* depending on specimen source may detect: *Adenovirus*, *CMV*, *enteroviruses*, *HSV*, *influenza*, *mumps*, *parainfluenza*, *RSV*, and *V-Z cell culture*

**Please Note:**

- **Multitest swab:** Up to **FOUR TESTS** can be run from **ONE** Aptima® Multitest Swab
- Specimens submitted from sources not noted as acceptable for the collection kit cannot be tested
- Do NOT include the cleaning swab in the collection vial; submit the Specimen Swab ONLY

| Women's Health Convenience Panels                                      |            |   |
|--|------------|---|
| SPECIMEN SOURCE<br>COLLECTION CONTAINER                                |            | VAGINAL<br>APTIMA® MULTITEST Swab (orange label)  |
| PANEL Name   | Order Code | Panel composed of:  |
| VAGINITIS PANEL (BV, CANDIDA 5 SPECIES, TV) NAA                        | LAB10638   | Total of 3 tests: <ul style="list-style-type: none"> <li>Candida Species Group, Candida Glabrata</li> <li>Trichomonas Vaginalis, NAA</li> <li>Bacterial vaginosis, NAA</li> </ul>   |
| VAGINITIS PLUS PANEL (BV, CANDIDA 5 SPECIES, TV, CT/GC) NAA            | LAB10639   | Total of 4 tests: <ul style="list-style-type: none"> <li>Candida Species Group, Candida Glabrata</li> <li>Trichomonas Vaginalis, NAA</li> <li>Bacterial vaginosis, NAA</li> <li>Chlamydia Trachomatis / Neisseria Gonorrhoeae RRNA</li> </ul> |
| CT/NG AND TRICHOMONAS NAA  | LAB10642   | Total of 2 tests: <ul style="list-style-type: none"> <li>Chlamydia Trachomatis / Neisseria Gonorrhoeae RRNA</li> <li>Trichomonas Vaginalis, NAA</li> </ul>  |
| Up to FOUR TESTS can be run from ONE Aptima® Multitest Swab Collection |            |   |

## Container Images

| APTIMA® MULTITEST | APTIMA® UNISEX | BD Affirm™ VPIII Swab | UTM® Viral Transport | eSwab®   |
|-------------------|----------------|-----------------------|----------------------|--|
|                   |                |                       |                      | <p><b>White cap:</b> Regular Flocked Swab<br/> <b>Purple cap:</b> Regular Flocked swab<br/> <b>Blue cap:</b> Flex Minitip Flocked swab<br/> <b>Green cap:</b> Minitip Flocked Swab</p> |