

URINE TOXICOLOGY SCREEN – FENTANYL NOW INCLUDED

March 22, 2022

Effective March 29, 2022, Pacific Diagnostic Laboratories (PDL) will be adding fentanyl to the Toxicology Screen, Urine LAB500.

TEST INFORMATION

Toxicology Screen (urine) effective March 29, 2022 will include the following drugs:

- Amphetamine Screen
- Barbiturate Screen
- Benzodiazepine Screen
- Cocaine Screen
- Opiate Screen
- Cannabinoid Screen
- Tricyclics Screen
- Urine PCP Screen
- Ecstasy Screen
- Fentanyl Screen

Will pricing of the Toxicology Screen (urine) change?

No. Pricing remains the same for the Toxicology Screen, urine.

Is the specimen collection a chain of custody (COD)?

No. PDL does not offer the collection as a chain of custody.

Is confirmation testing needed if a particular drug is detected on the screen?

The Toxicology Screen (urine) is a screening test performed using immunoassay test methodology. It is recommended that confirmation testing be done by mass spectrometry. Please note: PDL does not automatically reflex confirmation testing with this screen, an add-on request will be needed. Confirmation testing is sent out to our referral lab partner.

Test Detail:

Test Name	Order Code	Test Code	Order LOINC	CPT	CDM
Toxicology Screen, Urine	LAB500/UTOX	1230100479	69739-1	80307	30100484
Result name	Result Code	Result LOINC	Insurance Price		
AMPHETAMINE SCREEN, URINE [1534341]	1534341	19261-7	\$266.40		
BARBITURATES SCREEN, URINE [1534345]	1534345	19270-8			
BENZODIAZEPINE SCREEN, URINE [1558232]	1558232	14316-4			
COCAINE SCREEN, URINE [1811444]	1811444	19359-9			
OPIATES SCREEN, URINE [1555062]	1555062	19295-5			
CANNABINOIDS SCREEN, URINE [1558159]	1558159	18282-4			
TRICYCLICS SCREEN, URINE [1230001091]	1230001091	19312-8			
URINE PCP [1230017230]	1230017230	19659-2			
ECSTASY SCREEN, URINE [1230017229]	1230017229	19568-5			
FENTANYL SCREEN, URINE [1230018820]	1230018820	59673-4			
Specimen Requirements					
Specimen Source: Urine	Specimen Container: Sterile container	Transport Temperature: Refrigerate	TAT: 24 hours		

INTENDED USE:

Toxicology Screen (urine) LAB500 is designed for **PAIN MANAGEMENT or DRUGS OF ABUSE SCREENING**. It is not intended for workplace testing and does not comply with state regulatory workplace testing programs.

FOR MORE INFORMATION PLEASE CONTACT PDL CLIENT SERVICES AT 805-879-8100